


## MEMO 4D<sup>®</sup>

# Reshaping Mitral Repair

### ReChord: Chordal Guide System

Standardizes neochordae implantation




### Proven Annular Dynamics

Truly physiological three-dimensional motion to ensure natural behavior



### Innovation Starts at the Core

Unique Nitinol Cell Structure



### Stability and Support

Anterior Rigidity and Posterior flexibility to best support systolic function

### Exclusive Gradual Saddle Shaping

Designed to optimize the physiological geometry in enlarged annuli




42

### Unique Range of Sizes

Up to size 42 mm to treat large annuli

### Increase of the Anteroposterior Diameter

Helps accommodate excess leaflet tissue and reduces the risk of Systolic Anterior Motion (SAM)



MEMO 4D: DESIGNED TO BEST SUPPORT YOUR TECHNIQUE.

MEMO 4D's high-tech Nitinol core, thanks to its clinically proven physiological behavior, offers you a complete, long-lasting solution to treat MR. Now available with unique features to best support any repair.


# Information and accessories

## Product Ordering Information

Ordering Number	Size	A [mm]*	N° Loops
4DM-24	24 **	24	6
4DM-26	26	26	6
4DM-28	28	28	6
4DM-30	30	30	6
4DM-32	32	32	6
4DM-34	34	34	8
4DM-36	36	36	8
4DM-38	38	38	8
4DM-40	40	40	8
4DM-42	42	42	8

\* Inner metal core diameter

\*\* Made to order



## Accessories ordering information

Ordering Number	Name	Description
ICV0664	Uni Handle	Universal Bendable Handle
ICV1342	Extended Uni Handle	Universal Bendable Handle for MICS
ICV1357	Annuloplasty Ring Sizer Set	Complete Sizer Set (24 to 42 mm)
ICV1358	Annuloplasty Ring Accessory Tray	Empty Instrument Tray



Instrument Set



Uni Handle

Extended Uni Handle for MICS

MEMO 4D Sizer Set  
Now with more transparent heads for a better visibility



**LivaNova**

Health innovation that matters

[www.livanova.com](http://www.livanova.com)



Manufactured by:

**Sorin Group Italia Srl**

A wholly-owned subsidiary of LivaNova PLC  
Via Crescentino - 13040 Saluggia (VC) Italy  
Tel: +39 0161 487472 - Fax: +39 0161 487316  
[info.cardiacsurgery@livanova.com](mailto:info.cardiacsurgery@livanova.com)

MEMO 4D is FDA cleared.

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician. For indications, contraindications, precautions and warnings for each device, please refer to the Instructions For Use.

© 2018 LivaNova all rights reserved.